

Alco-Sensor IV with Memory

Alco-Sensor IV
with Memory



ALCOHOL TESTING INSTRUMENTS

- Evidential accuracy
- Automatic breath sampling
- Stores over 70 test results in memory
- Computer capable
- Printer capable
- DOT and law enforcement versions

ALCO
DRUG AND ALCOHOL TESTING PRODUCTS

1-800-227-9890

www.alcopro.com

Alco-Sensor IV

with Memory

Evidential Accuracy for Law Enforcement

The Alco-Sensor IV with Memory is approved by the U.S. Dept. of Transportation as an evidential breath testing device. The Alco-Sensor IV with Memory is widely used by law enforcement agencies as a PBT (Preliminary Breath Tester). Accuracy exceeds +/- .005 at .100 BrAC.

Dept. of Transportation Alcohol Testing

The Alco-Sensor IV with Memory connects to the DP1012 portable printer or the Slip printer to provide an economical and convenient DOT-compliant alcohol testing system. The Alco-Sensor IV with Memory assists the operator in complying with the DOT test protocol by counting down the 15-minute wait period after screening test results of .020 or greater and automatically starting a confirmation test.

Computer Capable

Law enforcement and corrections agencies use optional Total Recall software to facilitate uploading test results to a computer and generating reports from a database of test results.

Printer Options

DP1012 Printer



Power: Rechargeable batteries

Paper: Single ply or 3-ply, attach with tamper evident tape

Carrying Case: Pelican or Nylon case, also holds printer and AS IV with Memory

Slip Printer



Power: 12 volt AC

Paper: Prints directly onto Alcohol Test Form – no tamper evident tape

Carrying Case: No



SPECIFICATIONS

Sensor type: Electrochemical Oxidation (fuel cell)

Range of measurement: .000 - .400 BrAC

Accuracy and precision: U.S. DOT approved for evidential use

Breath sampling: Automatic with manual override

Response time: 5 – 10 seconds negative results, 10 – 30 seconds positive results

Power requirement: 9-volt battery (included)

Test mode: Direct

Calibration: Microprocessor controlled, using gas or wet bath simulator

Memory: More than 70 tests

Print capability: Yes, with DP1012 printer or Slip printer

Computer capability: Yes, Total Recall software or generic terminal emulation software

Temperature range: 0° - 40° C

Size: 3 1/2" x 5 1/4" x 1 1/2", 11 oz.

Warranty: One year, parts and labor

RELATED SUPPLIES

Item	Description	Price
312	Mouthpiece for Alco-Sensor IV	\$.24 @ 1-499 \$.23 @ 500-999 \$.22 @ 1,000-4,999 \$.20 @ 5,000-9,999 \$.17 @ 10,000+
324	Alcohol Gas Tank, 105 liter	\$165.00
360	Alcohol Gas Tank, 30 liter	\$139.00
328	Regulator	\$109.00
316	Using the Alco-Sensor IV DVD	\$50.00
347	Total Recall Software	\$525.00

Alco-Sensor IV with Memory

Item	Description	Price
206M	Alco-Sensor IV with Memory <i>includes hard-shell carrying case, 20 mouthpieces, calibration tool, battery, and manual.</i>	\$795.00
327-SLP	Slip Printer	\$680.00
213DP	Alco-Sensor IV with Memory and DP1012 Printer <i>includes Pelican carrying case, 100 mouthpieces, printer paper, printer ink ribbon, calibration tool, battery, battery charger, and manual.</i>	\$1,475.00
213DP-S	Alco-Sensor IV with Memory, DP1012 Printer, soft case <i>includes padded nylon carrying case, 30 liter tank, regulator, 100 mouthpieces, printer paper, printer ink ribbon, calibration tool, battery, battery charger, and manual.</i>	\$1,723.00

For Information & Orders Call: **1-800-227-9890**
 Order 24 Hours A Day At: **www.alcopro.com**
 Fax Orders To: **1-800-655-5834**

ALCO PRO
 DRUG AND ALCOHOL TESTING PRODUCTS